

Additional Mammographic Views: A Comprehensive Guide

Louise Miller, R.T.(R)(M)(ARRT), CRT(M), FSBI, FNCBC
Director of Education, Mammography Educators



MAMMOGRAPHY

© 2023 Mammography Educators

Additional Views Lexicon – Labeling Codes

- XCCL - Exaggerated craniocaudal lateral
- CV – Cleavage
- ML - Mediolateral
- LM - Lateromedial
- AT - Axillary tail
- TAN - Tangential



MAMMOGRAPHY

Additional Views Lexicon

- RL - Rolled lateral
- RM - Rolled medial
- RS - Rolled superior
- RI - Rolled inferior
- FB - Caudocranial



MAMMOGRAPHY

1

2

3


Additional Views Lexicon

- SIO - Superior lateral to inferior medial oblique
- LMO - Lateromedial oblique
- M - Magnification
- ID - Implant displaced
- No label: Spot Compression



MAMMOGRAPHY

4



MAMMOGRAPHY

Source: Hologic Selenia Unit

5

Labeling Codes (Lexicon)

The name is view (labeling code) is always preceded by identification of laterality

- i.e., LXCCL or RXCCL



MAMMOGRAPHY

6

The Most Commonly Used Additional Views

- XCCL
- CV
- LM/ML



MAMMOGRAPHY

7

Why do we do additional views?

- To show a specific component of the anatomy not seen on standard views
- To provide localization of an area of concern medial/lateral or superior/inferior to the nipple



MAMMOGRAPHY

8

OR. . .

- To show an area of concern in better detail
- To counteract superimposition of structures
- To triangulate a lesion



MAMMOGRAPHY

9

Commonly used additional views are done to show a specific component of the anatomy not seen on standard views.



10


Ask and Answer

- Which part of the breast do I want to visualize?
- In which projection?
- Which view will accomplish this?



11


Imaging the Breast in a Transverse or Axial Plane



12

XCCL - Exaggerated CC Lateral

Visualization of lateral breast tissue in a CC projection



13


Use of the XCCL

- Should be used on a baseline exam when lateral posterior breast tissue is missing on the CC view
- If glandular breast tissue on subsequent screening views is visualized back to the retromammary fat space on the MLO, an XCCL is not needed



14

XCCL required for baseline, but not on subsequent screenings



15


Use of the XCCL

- Should be performed on less than 10% of all patients
- Performed at 0 degrees angulation
- Patient's body should be at 45-degree angle to the IR
- Nipple should be pointing towards the upper corner of the image receptor



16


Muscle or NO muscle?
NO MUSCLE!




17




18




19




20




21




22




23




24



25



26



27

Axillary Tail



MAMMOGRAPHY

28

AT – Axillary Tail

- Angle tube to axillary tail
- Approximately 30 degrees



29

AT – Axillary Tail

It is never used to localize a lesion.



30

AT – Axillary Tail

- The AT View is used only for focal compression of the axillary tail
- Anterior to posterior orientation and compression



31

AT – Axillary Tail

It will not give you true lateral/medial or true superior/inferior orientation to the nipple



32


90-degree (True) Lateral

- LM - Lateromedial
- ML - Mediolateral



33

Superior or Inferior Orientation to the Nipple (LM or ML)



MAMMOGRAPHY

34

Use of the Lateral

- Shows effects of gravity on air fluid levels (Milk of Calcium)
- Used as a "tie breaker" view (to overcome superimposition of structure)
- Visualizes the breast in the sagittal plane (demonstrates an area of concern superior or inferior to the nipple)



35

Why do the LM?

- When you did the MLO, you showed the lateral breast in better detail; The LM shows the medial breast in better detail
- The LM takes advantage of the lateral mobile border of the breast and thus facilitates positioning



36

Why do the LM?

- The posterior medial breast is hardest part of the breast to image and the area most often missed on the MLO
- If done properly, by off-setting the IR into the contralateral breast, you will be able to go deeper against the chest wall



37


Why do the LM?

- There is no issue of the contralateral breast impeding the path of the compression paddle




38

Improperly positioned LM with breasts separated, so the middle of the IR is centered on midsternal line. This excludes deep medial breast tissue on the side you are imaging.



39

Properly positioned LM with breasts separated so the top edge of the IR is centered on midsternal line and the width of the IR pressing against the contralateral breast.



40

ML – Mediolateral

The opposite breast must be pulled back to allow the compression paddle to pass and may therefore eliminate visualization of deep medial breast tissue.



41

Additional Views to Facilitate Imaging of Difficult Body Habitus

- LMO
- FB



42

LMO – Lateromedial Oblique FB – Caudocranial (From Below)

- Used when a standard MLO or CC is impossible
- Kyphotic patients
- Males with prominent pectoral muscles



43

LMO – Lateromedial Oblique


- Degree of angulation should be 90 degrees plus or minus the standard degree of angulation.
- Average 45-degree RMLO would be 135 degrees for RLMO



44

MLO

LMO



45

FB - Caudocranial

- Machine is turned 180 degrees opposite the CC
- Direction of the beam is Caudal to Cranial




MAMMOGRAPHY

46

CC

FB



MAMMOGRAPHY

47

Additional Views for Clarification of Areas of Concern

- TAN
- Spot compression
- Spot compression with MAG
- Rolled views



MAMMOGRAPHY

48


TAN – Tangential View

- To prove the existence of dermal calcifications
- Enhanced visualization of palpable masses that may otherwise be superimposed on glandular breast tissue




MAMMOGRAPHY

49



MAMMOGRAPHY

50




MAMMOGRAPHY

51



MAMMOGRAPHY

52



MAMMOGRAPHY

53

Localization for Verification of Skin Calcifications

- Decreased with use of DBT
- Set up as the same as a needle localization
- Determine which quadrant the calcifications are located




MAMMOGRAPHY

54



Localization for Verification of Skin Calcifications

- Use biopsy paddle
- Select direction of approach so that the window of biopsy paddle is closest to the area in question





55

Tangential Views for the Verification of Skin Calcifications

56

Tangential Views for the Verification of Skin Calcifications

57






Image Courtesy of Robyn Hadley, R.T.(R)(M)



58

Tangential Views for the Verification of Skin Calcifications




59

Spot Compression with/without Magnification



60

All Would be Labeled LTAN






61

Spot Compression Paddles

62

Square vs Round Paddle

63


Spot Compression with Magnification



MAMMOGRAPHY

64

Square vs Round Paddle





MAMMOGRAPHY

Image Courtesy of Robyn Hadley, R.T.(R)(M)


65

Spot/Mag Measurements



MAMMOGRAPHY

66





MAMMOGRAPHY

Image Courtesy of Robyn Hadley, R.T.(R)(M)

67



Spot/Mag Measurements



MAMMOGRAPHY

68

Spot/Mag Measurements



MAMMOGRAPHY

69

Remember...

- You must simulate compression when making measurements on the breast
- Mark the center of the target area with a surgical marker so you can make appropriate corrections on subsequent images, if needed



71

Remember...

- Roll views** – Used to overcome superimposition of structures by changing the orientation of the beam to the breast
- Lateral views** – Used to overcome superimposition of structures by changing the orientation of the breast to the beam



72



MAMMOGRAPHY

Image Courtesy of Robyn Hadley, R.T.(R)(M)

70

Rolled Views


- RCCRM – Right CC superior breast rolled medial
- RCCRL – Right CC superior breast rolled lateral



MAMMOGRAPHY

73


LCCRM – Left CC Rolled Medial



MAMMOGRAPHY

74


Lateral Used to Overcome Superimposition



MAMMOGRAPHY

75

Lateral Used for Lesion Location



MAMMOGRAPHY


76

Roll View to Overcome Superimposition

1. Something seen on CC, but not on MLO - is it real?



2. Do LCCRM:



Examination is complete - this was superimposition and not a "real" mass.



MAMMOGRAPHY

77

Roll View for Lesion Location

Example A:



Example B:



MAMMOGRAPHY

78

Imaging of Augmented Breasts

- CCs views of each breast with implants in place
- MLO views of each breast with implants in place
- CCID views of each breast with implant displaced
- MLOID views of each breast with implant displaced



MAMMOGRAPHY

79


Imaging of Augmented Breasts: Full Implant Views

- Should be done with only enough compression to immobilize the breast to prevent motion unsharpness
- Curved paddle can be used (if available)
- Appropriate technique (usually manual) should be used



MAMMOGRAPHY

80



MAMMOGRAPHY

81

Imaging of Augmented Breasts: ID Views

- ID views (depending on implant mobility) can be performed with taut compression
- Half paddle can be used for patients with small amount of natural breast tissue
- Appropriate technique (used for patients without implants) should be used (implants)
- Patient can be positioned from behind (with tech standing and/or patient seated)



MAMMOGRAPHY

82


Half Paddle



MAMMOGRAPHY

83


Imaging of Augmented Breasts



MAMMOGRAPHY

84

ID – Implant Displaced



MAMMOGRAPHY

Source: lebradictionary.com

85


Working from Behind for CCID Views



MAMMOGRAPHY

86

Working from Behind for MLOID Views



MAMMOGRAPHY

87

Conclusion

- Additional views are helpful in identifying true location of areas of concern
- They are used for diagnostic work-ups
- Can provide valuable information to aid in diagnosis of breast cancer



MAMMOGRAPHY

88

References

Miller, L. C., RTRM. (2015). *Mammography Positioning Guidebook*. CA: Mammography Educators.

Miller, L. C., RTRM (ARRT), Lohmann, T. L., B.S. & T. (ARRT). (2020). *Image Quality & Positioning Problem-Solving for Breast Imaging: Meeting EQUIP Standards* (1st ed.). San Diego: California Mammography Educators.

Long, S. M., Miller, L. C., Bosco, M. A., & Martin, L. L. (2010). *Handbook of Mammography* (5th ed.). Edmonton: Mammography Consulting Services.

Derenberger, Dawn, Hadley, Robyn. "Most Commonly Used Additional Views, Part 3: Defining Structures and Clarifying Presence of Abnormalities." *SBI News* Issue 3, 2021. <https://www.sbi-news.com/most-commonly-used-additional-views-part-3-defining-structures-and-clarifying-presence-of-abnormalities>



MAMMOGRAPHY

89

For questions or more

www.information.com

619-663-8269

info@mammographyeducators.com

Services we offer, include:

- Onsite Positioning Training
- Assistance with Accreditation & Inspection Processes
- Live Webinars
- Customized Continuing Education Programs



MAMMOGRAPHY

90