

MRI-Guided Breast Biopsy

Terry Lehmann, B.S., R.T.(R)(M)
Breast Imaging Consultant, Mammography Educators



MAMMOGRAPHY
EDUCATORS

© 2025 Mammography Educators

1

MRI of the Breast

- MRI is an important tool for the detection of breast cancer and assessment of silicone implant integrity
- The use of MRI with contrast agents for the detection of breast cancer was first reported in 1986



MAMMOGRAPHY
EDUCATORS

Heywang 1986

2

Current Indications for Breast MRI

- Axillary metastasis of unknown origin
- High risk patients (screening MRI)
- Extent of disease (with new diagnosis)
- Aid to pre-surgical treatment planning
- Treatment response/follow-up
- History of chest radiation



MAMMOGRAPHY
EDUCATORS

3

Current Indications for Breast MRI

- Differentiating postoperative changes from recurrence
- Reduced mammographic sensitivity
- Problem solving for selective inconclusive clinical or mammographic/ultrasound findings
- Guidance for biopsy



MAMMOGRAPHY
EDUCATORS

4

Breast MRI Technique

- MRI creates images of the breast by measuring changes in the movement of protons in fat and water with the application of changing magnetic fields
- By utilizing the differences in tissue relaxation characteristics, an image is acquired by processing the signal changes that occur following application of pulses of energy



5

MRI for Breast Cancer Detection

- Requires administration of the contrast agent gadolinium
- The use of contrast is based on the concept of tumor angiogenesis or neovascularity



6

Breast MRI Technique

Tumor-associated blood vessels have increased permeability, which leads to prompt uptake and release of gadolinium within the first 1-2 minutes after administration, leading to a pattern of rapid enhancement and washout on MRI.



7

MR Breast Lexicon

- Focus/foci (tiny spot <5mm) no other descriptor
- Mass (3D space-occupying lesion)
 - Shape (round, oval, lobulated, irregular)
 - Margin (smooth, irregular, spiculated)
 - Mass enhancement (homogeneous, heterogeneous, rim, non-enhancing septations, enhancing septations)



8

MR Breast Lexicon

- Non-mass enhancement
 - Focal, linear, ductal, segmental, regional, multiple regions, diffuse
 - Homogeneous, heterogeneous, stippled/punctate, clumped NMLE
- Symmetric or asymmetric
- Other: nipple retraction, edema, lymphadenopathy, pectoralis muscle or chest wall invasion, skin thickening/invasion
- Kinetic curves



MAMMOGRAPHY
EDUCATORS

9

MRI of the Breast

Evaluates:

- Morphology
- Kinetics/Angiogenesis



MAMMOGRAPHY
EDUCATORS


10

Kinetic Curves

Washout

Plateau


Persistent



MAMMOGRAPHY
EDUCATORS

11

MRI Findings Flow Chart



MAMMOGRAPHY
EDUCATORS

12

Breast MRI

- Correlate with mammogram/ultrasound (US)
- Second look mammogram/US after MRI imaging/finding
- May be able to convert biopsy to US guidance



MAMMOGRAPHY
EDUCATORS

13

MRI Breast Biopsy – the “Why”


- Performed if lesion cannot be identified by other modalities
- 14g core can be done
- Vacuum assisted core is gold standard
- Dedicated coil for biopsy
- Varies from standard bilateral coil used
- Typically, far fewer coils
- Unilateral imaging only



MAMMOGRAPHY
EDUCATORS

14

“Our” Standard Coil – 16 Channel



MAMMOGRAPHY
EDUCATORS

Courtesy Philips

15

“Our” Biopsy Coil – 7 Channel



MAMMOGRAPHY
EDUCATORS

Courtesy Philips

16

Pitfalls

- Some systems are limited to lateral approach only
- Breast size
- Posterior lesions
- IMF lesions
- High axillary lesions



MAMMOGRAPHY
EDUCATORS

17




MAMMOGRAPHY
EDUCATORS

18



MAMMOGRAPHY
EDUCATORS

19



MAMMOGRAPHY
EDUCATORS

20



21

Patient JMF

- Hx prior Stage IIA IDCA 2005
- Triple neg
- Stage
- BRCA1+
- Treated with BCS, chemo & XRT
- Annual MRI added to surveillance



22

“The Why”

Finding L1 CAD Summary


Location
L: LO, 5 o'clock, middle

Cor Ax-CC Sag-ML

L Med Lat

Distance (cm) edge to:
Nipple: 3.9 (radial 1.9)

Size
Diameters: 1.1 x 0.42 cm



23


MRI Breast Biopsy – The “How”

- Position patient in biopsy cradle
- Skin is cleaned *before* patient is positioned in the biopsy coil
- It will be cleaned a 2nd time, before giving anesthetic



24

- Prone, feet first
- Comfortable headrest, pillow



MAMMOGRAPHY
EDUCATORS

25

MRI Breast Biopsy Position

- Breast is supported in the biopsy cradle/coil by 2 plates
- Grid on aspect/side of breast accessed for biopsy stabilizes breast
- Opposite side of breast supported by solid plate



MAMMOGRAPHY
EDUCATORS

26


Grid and Solid Compression Plates



MAMMOGRAPHY
EDUCATORS

27


Solid Plate "Unlocked"



MAMMOGRAPHY
EDUCATORS

28


Solid Plate “Locked”



MAMMOGRAPHY
EDUCATORS

29

Grid Support Plate “Unlocked”



MAMMOGRAPHY
EDUCATORS

30

Introducer Kit and Tissue Marker “Clip”



MAMMOGRAPHY
EDUCATORS

31


Fiducial



MAMMOGRAPHY
EDUCATORS

32

Fiducial Placed in Grid



33

MRI Breast Biopsy – The “How”

- Perform “scout” scan for parameters
- Give GAD injection
- Perform dynamic imaging
- “Send” images to CAD program
- Localize skin fiducial & lesion for targeting data, using CAD program



34


MRI CAD Program

- Provides lesion identification, analysis and kinetics from screening or work-up study
- Used with biopsy procedures for targeting lesion or ROI to be sampled




35

Patient Positioned in BX Coil/Cradle




36

Patient Positioned in BX Coil/Cradle




37

First Post-contrast Image



38

Skin Fiducial



39

CAD Targeting Data

Localizer 1: Lateral view of right breast

Localizer 1


Move relative to the reference point:		
Distance	Direction	Toward
3 blk* (5.4 cm)	Anterior	nipple
2 blk* (3.1 cm)	Superior	Head
5.8 cm	Depth	

* Assumes Reference is in center of grid block
 Insertion: Lateral Right
 Method: SURCOS BIOPSY



40


CAD Targeting Data



MAMMOGRAPHY
EDUCATORS

41

CAD Targeting Data



MAMMOGRAPHY
EDUCATORS

42

MRI Breast Biopsy – The “How”

- Patient's skin is cleaned again, through the grid section identified by targeting
- Universal time outperformed
- Patient given superficial and deep lidocaine
- Skin nick made with scalpel blade
- Plastic introducer stylet and trocar placed in breast (with depth set on hub)



MAMMOGRAPHY
EDUCATORS

43




MAMMOGRAPHY
EDUCATORS


44




45



46




47




48

Placement of Needle Guide



MAMMOGRAPHY
EDUCATORS


49



MAMMOGRAPHY
EDUCATORS

50


Needle Guide Used to Support Introducer and Trocar



MAMMOGRAPHY
EDUCATORS

51


Obturator Placed Through Needle Guide



MAMMOGRAPHY
EDUCATORS

52

Final Position Flush Against Needle Guide



MAMMOGRAPHY
EDUCATORS

53

MRI Breast Biopsy – The “How”

- Plastic obturator replaces trocar
- Needle confirmation scan done
- To ascertain that needle/obturator tip is at desired location relative to target area
- If accuracy confirmed, proceed with tissue sampling



MAMMOGRAPHY
EDUCATORS

54

Needle Confirmation Scan



MAMMOGRAPHY
EDUCATORS

55

MRI Breast Biopsy – The “How”


- Confirmation scan is reviewed for accuracy and need for any adjustments
- Plastic obturator is removed
- Tissue sampling can begin
- VAC assisted probe inserted through needle guide to targeted lesion



MAMMOGRAPHY
EDUCATORS

56


Biopsy Probe



MAMMOGRAPHY
EDUCATORS

57

Probe Replaces Obturator



MAMMOGRAPHY
EDUCATORS

58

Biopsy Samples Taken


- Probe rotated around clockface
- Quadrant or hemisphere or “clock” identified by relative position of ROI/lesion to probe
- Vacuum assisted sampling
- Lavage after sampling



MAMMOGRAPHY
EDUCATORS

59

9 g Vacuum Assisted System



Stays outside of room-
Console NOT MRI compatible



MAMMOGRAPHY
EDUCATORS

60

Tissue Samples Retrieved



61


MRI Breast Biopsy – The “How”

- Marker clip placement
- Inserted through introducer, similar to stereotactic biopsy
- Post-clip placement scan may be done if desired



62

Post-clip Scan



63


MRI Breast Biopsy – The “How”

After hemostasis is achieved, wound dressing placed and post-biopsy mammogram done for clip placement evaluation.



64

Post-biopsy Mammogram



MAMMOGRAPHY
EDUCATORS

65

Pathology Results

- DCIS intermediate grade
- Concordant with Imaging findings
- Patient elected bilateral mastectomy
- SLNB negative
- Patient went on to 2-stage reconstruction



MAMMOGRAPHY
EDUCATORS

66

MRI Breast Biopsy – The “How”

- Position patient (after cleaning skin)
- Perform “scout” scan for imaging parameters
- Give GAD injection
- Perform dynamic imaging
- Localize skin fiducial and lesion for targeting data with CAD program
- Clean skin a second time
- Universal time-out pause



MAMMOGRAPHY
EDUCATORS

67

MRI Breast Biopsy – The “How”

- Administer superficial & deep lido
- Make skin nick
- Insert needle guide into grid
- Insert introducer and advance trocar to target
- Replace trocar with obturator
- Perform needle confirmation scan



MAMMOGRAPHY
EDUCATORS

68

MRI Breast Biopsy – The “How”

- Remove obturator and perform tissue sampling
- Place marker clip
- Do post biopsy mammogram
- Follow-up as directed by pathology results



MAMMOGRAPHY
EDUCATORS

69

Summary

- MRI is an important tool for the detection of breast cancer and assessment of silicone implant integrity
- Images are analyzed for morphology and initial/delayed contrast enhancement
- Sensitivity is high, but specificity is lower
- Provides biopsy guidance for MRI-only imaging findings



MAMMOGRAPHY
EDUCATORS

70

References

- Chernoff Daniel, MD PhD, *Principles of magnetic resonance imaging*. UpToDate. Jan 2020.
- Heywang SH, Hahn D, Schmidt H, et al. *MR Imaging of the breast using gadolinium*. J Computer Assisted Tomography 1986; 10:199
- Lehman CD, Gatsonis C, Kuhl CK, et al. *MR evaluation of the contralateral breast in women with recently diagnosed breast cancer*. NEJM 2007; 356:1295
- Slanetz, Patricia J, MD,MPH,FACR. *MR of the breast and emerging technologies*. UpToDate. Aug 2020.



MAMMOGRAPHY
EDUCATORS

71

Thank You!

Services we offer, include:

- Onsite Positioning Training
- Assistance with Accreditation & Inspection
- Live Webinars and Conferences
- On-Demand Continuing Education

For questions or more information:

619-663-8269

mammographyeducators.com

info@mammographyeducators.com



MAMMOGRAPHY
EDUCATORS

72